

TrackCore Operating Room, the industry leading tissue and implant management software, integrates with intelligent RFID-enabled enclosures to automate your chain of custody and inventory management without human intervention.

AUTOMATE CHAIN OF
CUSTODY TRACKING

IMPROVE REGULATORY
COMPLIANCE

ABSOLUTE
ACCOUNTABILITY

INVENTORY
OPTIMIZATION



HOW RFID TECHNOLOGY WORKS

ONE

Place the internet-enabled RFID enclosure at the point of use.

TWO

Pre-authorized users are able to access the enclosure with RFID-enabled badges.

THREE

Users can then add or remove RFID-tagged inventory from the enclosure.

FOUR

The user, items added or removed, date/time, and temperature are all automatically tracked.



TC-179 | ULT Under Counter Freezer | 3.6ft³ Intelligent Freezer

DETAILS AND DIMENSIONS

COMPONENT	SPECIFICATION
External Dimensions	32.7"H x 35.4"W x 25"D Place 3" from back wall Allow 21" for door clearance
Internal Dimensions	25.2"H x 18.9"W x 13"D
Volume	3.6 ft ³
Weight	265 lbs
Temperature Range	- 86° C to - 50° C
Security Device	Electronic Door Lock
Access Control	LF pass reader
Refrigerant	CFC/HCFC-free Refrigerant Mixtures/Non-flammable
Cooling Method	2-stage cascade system, CO2 & LN2 Backup Available



COMMUNICATION	SPECIFICATION	POWER	SPECIFICATION
RFID Frequencies	902 to 928 MHz	Electrical Power	120 VAC, 60 Hz, 13 A
Communications	Ethernet, optional cellular	Power Consumption	7.2 kWh/day
		RFID Battery Backup	24 VDC, 1.3 Ah sealed lead acid

W trackcoreinc.com
 P 616.632.2222
 E sales@trackcoreinc.com

CONFIDENTIAL: This document is intended only for the original recipient. The information contained herein is confidential and may constitute inside or non-public information under international, federal, or state securities laws. Unauthorized forwarding, printing, copying, distribution or use of such information is strictly prohibited and may be unlawful. Please contact TrackCore, Inc. at +1.616.632.2222.

