

TrackCore Operating Room, the industry leading tissue and implant management software, integrates with intelligent RFID-enabled enclosures to automate your chain of custody and inventory management without human intervention.

AUTOMATE CHAIN OF  
CUSTODY TRACKING

IMPROVE REGULATORY  
COMPLIANCE

ABSOLUTE  
ACCOUNTABILITY

INVENTORY  
OPTIMIZATION



### HOW RFID TECHNOLOGY WORKS

ONE

Place the internet-enabled RFID enclosure at the point of use.

TWO

Pre-authorized users are able to access the enclosure with RFID-enabled badges.

THREE

Users can then add or remove RFID-tagged inventory from the enclosure.

FOUR


The user, items added or removed, date/time, and temperature are all automatically tracked.



# TC-196 | Undercounter RFID Freezer | 4.2ft<sup>3</sup> Intelligent Freezer

## DETAILS AND DIMENSIONS

COMPONENT	SPECIFICATION
External Dimensions	34"H x 23.75"W x 24"D
Internal Dimensions	24"H x 19"W x 17.5"D
Volume	4.2 ft <sup>3</sup>
Weight	140 lbs
Drawer Capacity	10.0 lbs. per drawer
Temperature Range	- 25° C to - 15° C
Security Device	Electronic Door Lock
Access Control	LF Pass Reader
Refrigerant	HFC-Free R600a
Cooling Method	In-shelf tubing



COMMUNICATION	SPECIFICATION	POWER	SPECIFICATION
RFID Frequencies	902 to 928 MHz	Electrical Power	115 VAC, 60 Hz, 2 A
Communications	Ethernet, optional cellular	Power Consumption	1.2 kWh/day
		RFID Battery Backup	24VDC, 1.3 Ah sealed lead acid

W [trackcoreinc.com](http://trackcoreinc.com)
P 616.632.2222
 E [sales@trackcoreinc.com](mailto:sales@trackcoreinc.com)

CONFIDENTIAL: This document is intended only for the original recipient. The information contained herein is confidential and may constitute inside or non-public information under international, federal, or state securities laws. Unauthorized forwarding, printing, copying, distribution or use of such information is strictly prohibited and may be unlawful. Please contact TrackCore, Inc. at +1.616.632.2222.

